


how to choose...



	Class I	Class II		Class III	
		Moderate	Severe	Very severe	
3 y.o.			EF Kid 		
6 y.o.		EF Start 		EF Classe III Petit 	
8 y.o.	EF T Slim  EF Profil 	EF T Slim 	EF Classe II Slim 	EF Classe II two steps 	EF Classe III Petit  EF Classe III Standard 
10 y.o.	EF Profil  EF Guide 	EF T Slim Long 	EF Classe II Standard  Deep labiomental fold	 Treatment of overjet greater than 10 mm	EF Classe III Standard 
12 y.o.		EF Guide 			
and +					

Warning: This selection is for information only. Arch size, dental age and treatment goals vary in each patient.

EF T Slim




CL I

CL II
moderate

8 - 10
y.o.

3 to 3 indents



Physiological rest position

- Universal 3 to 3 maxillary and mandibular indents to align the midlines.

Available in 5 colors:



D590 300R

D590 300B

D590 300V

D590 300J

D590 300

Case treated by EF T Slim

Martin 10 y.o., before EF T Slim

- Class II 2 moderate



Martin, after 12 months EF T Slim



EF T Slim Long




CL I

CL II
moderate

10 - 12
y.o.

3 to 3 indents



Physiological rest position

- Universal 3 to 3 maxillary and mandibular indents to align the midlines.
- Distal extension to include 6s and 7s.

Available in 5 colors :



D590 300LR

D590 300LB

D590 300LV

D590 300LJ

D590 300L

Case treated with EF T Slim Long

Aurélien 10 y.o., before EF T Slim Long

- Skeletal Class I • Tongue dysfunction when swallowing • Family history (mother had maxillofacial surgery and father was treated for impacted maxillary canines)



Aurélien, after 12 months EF T Slim Long

